

# Microfluidic Immunohistochemistry for Rapid Detection of Diagnostic and Predictive Biomarkers on Non-Small Cell Lung Carcinoma

Procopio Maria-Giuseppina<sup>1</sup>, Brajkovic Saska<sup>2</sup>, Leblond Anne-Laure<sup>1</sup>, Dupouy Diego Gabriel<sup>2</sup> and Soltermann Alex<sup>1</sup>

1. Institute of Pathology and Molecular Pathology, University Hospital Zurich, CH-8091 Zurich, Switzerland | 2. Lunaphore Technologies SA, EPFL Innovation Park-Building C, CH-1015 Lausanne, Switzerland

## Background

Lung cancer is the most lethal tumor world-wide. Immunohistochemistry (IHC) is an essential tool in oncologic pathology. Considering the increasing cancer incidence, a concomitantly accurate and time-efficient IHC is clinically highly relevant. The automated microfluidic staining device enables for a precisely controlled immune-reaction to take place in an extremely short incubation time; typically, a complete microfluidic IHC lasts less than 12 minutes on frozen tissue section.

## Objective

We aimed to develop and validate a microfluidics IHC for fast and accurate biomarkers detection on formalin fixed paraffin-embedded (FFPE) non-small cell lung carcinomas (NSCLC) - the most common lung cancer type. Here, we focus on the first part of the study aiming at optimizing the microfluidic protocols.


## Methods

Each protocol step - from epitope retrieval to counterstaining - was automatized and optimized on the microfluidic staining device (**Figure 1**) for all biomarkers on representative FFPE NSCLC specimens to reach high-quality staining equivalent to a routinely used Ventana BenchMark ULTRA automated stainer.

## Results

The optimized microfluidic IHC protocol achieved an analytical performance comparable to standard staining on NSCLC for (i) differential diagnosis (**Figure 2**) into adenocarcinoma (TTF1/CK7) or squamous cell carcinoma (p40/CK5-6) and (ii) prediction of immuno-therapy response (PD-L1) (**Figure 3**). Concomitantly, the total process time was shortened down to twenty minutes - including the epitope retrieval pre-staining step, with a time reduction up to nine folds (**Table 1**).

**FIGURE 1.** Microfluidic device schematic representation




**TABLE 1.** Optimized microfluidic IHC protocols in comparison with standard assays on FFPE NSCLC


Marker	p40		TTF1		CK5-6		CK7		PD-L1	
	Abl Clone	BC28	SP141	SP141	D5/16B4	D5/16B4	SP52	SP52	E1L3N	IHC411
Abl Dilution	1:100	1:50	1:50	RTU	RTU	RTU	RTU	RTU	1:100	1:50
Staining Device	STD	MTP	STD	MTP	STD	MTP	STD	MTP	STD	MTP
HIER	24	5	56	5	36	5	4*	10	56	10
H <sub>2</sub> O <sub>2</sub>	4	--	4	--	4	--	4	--	4	--
Abi	32	4	12	4	32	1	16	1	32	4
Linker	8	--	8	--	8	--	8	--	8	--
AbII	8	4	8	4	8	1	8	1	8	4
DAB	8	2	8	2	8	2	8	2	8	2
CuSO <sub>4</sub>	4	**	4	--	4	--	4	--	4	**
HTX	8	0.25	8	0.25	8	0.25	8	0.25	8	0.25
<b>Total Time (min)</b>	<b>173</b>	<b>26</b>	<b>184</b>	<b>26</b>	<b>188</b>	<b>20</b>	<b>124</b>	<b>25</b>	<b>220</b>	<b>31</b>

Optimized microfluidic IHC protocols (MTP) for p40, TTF1, CK5-6, and CK7 on FFPE NSCLC. Incubation time for all protocol step and primary antibody dilutions are shown as compared to routinely used standard (STD) protocols. RTU: ready to use; HIER: heat-induced epitope retrieval; H<sub>2</sub>O<sub>2</sub>: oxygen peroxidase; Abi/II: primary/secondary antibody; DAB: 3,3'-diaminobenzidine; CuSO<sub>4</sub>: copper sulfate; HTX: hematoxylin. \*PIER: proteolytic-induced epitope retrieval \*\*CuSO<sub>4</sub> post-staining manual incubation for 30sec.

**FIGURE 2.** Analytical performance comparison for p40, TTF1, CK5-6 and CK7 IHC staining by standard versus microfluidic assay on FFPE NSCLC tissues



**FIGURE 3.** Analytical performance comparison for PD-L1 IHC staining by standard versus microfluidic assay on FFPE tonsil and NSCLC tissues



## Conclusions and Future Directions

The microfluidic IHC resulted in fast automated high-quality staining for all assessed markers on NSCLC. Subsequent clinical validation on a large cohort will provide a diagnostic tool for biomarkers detection in a turnaround time far beyond the existing automated stainers.